

# Mini Spin Column Systems



## CUSTOMISED SOLUTIONS FOR

FILTRATION

ULTRA FILTRATION

DNA | RNA PURIFICATION



**Mini Spin Column**  
round bottom  
without lid

- Developed for the use in purification, isolation and separation of biomolecules, particularly DNA/RNA
- Filter tubes can be delivered with various membrane filters (e.g. CA|RC|GF|PVDF|PE|PTFE|Nylon) depending on customer specific applications\*
- Available with 1.5 ml or 2.0 ml receiver tubes



**Mini Spin Column**  
straight bottom  
with lid



\*For further information or product specific inquiries contact us [info@ahn-bio.de](mailto:info@ahn-bio.de)

### Spin Column Single Parts without filter

Spin Column Single Parts without filter	Cat. No.	Pack Type	Sterilized	Sales Unit
0.8 ml, Spin Column Filter Tube, clear	3-330-50-0	Bag		1x1000
0.8 ml, Spin Column Filter Tube, clear, with filtering	3-331-50-0	Bag		1x1000
0.8 ml, Spin Column Filter Tube, clear, with cap	3-336-50-0	Bag		1x1000
1.5 ml, Receiver Tube, clear	3-200-50-0	Bag		1x1000
	3-200-80-0	Bag		10x1000
2.0 ml, Receiver Tube, clear	3-202-50-0	Bag		1x1000
	3-202-80-0	Bag		10x1000