

# **COMPATIBILITY OF**

AHN myTip LT Low Retention Tips

# AHN myTip Low Retention Tips Compatibility with AHN Single Channel Pipettes

Tips Compatibility  
AHN Single Channel Pipettes

Pipettor		AHN																														
Brand		pipet4u <sup>®</sup> pro										pipet4u <sup>®</sup> Revolution																				
Model		variable/single channel					fix/single channel					variable/ electronic single channel																				
Volume Range / Cat. No.		0.1-2.5 µl / 8-100-00-9	0.5-10 µl / 8-101-00-9	2-20 µl / 8-102-00-9	5-50 µl / 8-103-00-9	10-100 µl / 8-104-00-9	20-200 µl / 8-105-00-9	50-200 µl / 8-105-01-9	100-1000 µl / 8-106-00-9	200-1000 µl / 8-106-01-9	500-5000 µl / 8-107-00-9	2000-10000 µl / 8-108-00-9	5 µl / 8-001-00-9	10 µl / 8-002-00-9	20 µl / 8-003-00-9	25 µl / 8-004-00-9	50 µl / 8-007-00-9	100 µl / 8-008-00-9	200 µl / 8-009-00-9	250 µl / 8-010-00-9	500 µl / 8-011-00-9	1000 µl / 8-012-00-9	2000 µl / 8-013-00-9	2500 µl / 8-014-00-9	5000 µl / 8-015-00-9	10000 µl / 8-016-00-9	0.2-10 µl / 8-500-00-0	5-100 µl / 8-501-00-0	10-200 µl / 8-502-00-0	50-1000 µl / 8-503-00-0	100-5000 µl / 8-505-00-0	
		AHN myTip LT Tips																														
Volume	Cat.No.																															
0.1-10 µl	4-000-81-0 4-000-C5-0		■	■	■																							■				
0.1-10 µl XL	4-060-81-0 4-060-C5-0		■	■	■																							■				
0.5-10/20 µl	4-050-81-0 4-050-C5-0			■																												
1-200 µl	4-110-81-0 4-110-C5-0				■	■	■	■				■	■	■	■	■	■											■	■			
100-1000 µl long	4-200-80-0 4-200-C5-0								■	■										■	■										■	■



# AHN myTip Low Retention Tips Compatibility with AHN Multi Channel Pipettes

Pipettor		AHN											
Brand		pipet4u <sup>®</sup> pro				pipet4u <sup>®</sup> Revolution							
Model		variable/multi channel				variable/multi channel							
Volume Range / Cat. No.		8 Channel		12 Channel		8 Channel		12 Channel					
		0.5-10 µl / 8-20100-9	5-50 µl / 8-202-00-9	30-300 µl / 8-203-00-9	0.5-10 µl / 8-401-00-9	5-50 µl / 8-402-00-9	30-300 µl / 8-403-00-9	0.2-10 µl / 8-520-00-0	2.5-50 µl / 8-521-00-0	15-300 µl / 8-522-00-0	0.2-10 µl / 8-530-00-0	2.5-50 µl / 8-531-00-0	15-300 µl / 8-532-00-0
AHN myTip LT Tips													
Volume	Cat.No.	8 Channel		12 Channel		8 Channel		12 Channel		8 Channel		12 Channel	
0.1-10 µl	4-000-81-0 4-000-C5-0	■	■	■	■	■	■	■	■	■	■	■	■
0.1-10 µl XL	4-060-81-0 4-060-C5-0	■	■	■	■	■	■	■	■	■	■	■	■
0.5-10/20 µl	4-050-81-0 4-050-C5-0	■	■	■	■	■	■	■	■	■	■	■	■
1-200 µl	4-110-81-0 4-110-C5-0	■	■	■	■	■	■	■	■	■	■	■	■

**Tips Compatibility**  
AHN Multi Channel Pipettes





Uthleber Weg 14  
D-99734 Nordhausen  
Germany

Tel.: +49(0)3631/46594-04  
Fax: +49(0)3631/46594-10

E-Mail: [info@ahn-bio.de](mailto:info@ahn-bio.de)  
Web: [www.ahn-bio.de](http://www.ahn-bio.de)