

# KOMPATIBILITÄT

AHN myTip LT Low Retention Pipettenspitzen



# AHN myTip Low Retention Pipettenspitzen Kompatibilität mit AHN Mehrkanal-Pipetten

Pipette		AHN											
Marke		pipet4u <sup>®</sup> pro				pipet4u <sup>®</sup> Revolution							
Modell		variabel/Mehrkanal				variabel/Mehrkanal							
		8-Kanal		12-Kanal		8-Kanal		12-Kanal					
Volumenbereich / Kat.-Nr.		0.5-10 µl / 8-201-20-9	5-50 µl / 8-202-20-9	30-300 µl / 8-203-20-9	0.5-10 µl / 8-401-20-9	5-50 µl / 8-402-20-9	30-300 µl / 8-403-20-9	0.2-10 µl / 8-520-00-0	2.5-50 µl / 8-521-00-0	15-300 µl / 8-522-00-0	0.2-10 µl / 8-530-00-0	2.5-50 µl / 8-531-00-0	15-300 µl / 8-532-00-0
AHN myTip LT Pipettenspitzen													
Volumen	Kat.-Nr.												
0.1-10 µl	4-000-81-0	■			■			■				■	
	4-000-C5-0	■			■			■				■	
0.1-10 µl XL	4-060-81-0	■			■			■				■	
	4-060-C5-0	■			■			■				■	
0.5-10/20 µl	4-050-81-0	■			■			■				■	
	4-050-C5-0	■			■			■				■	
1-200 µl	4-110-81-0	■	■		■	■		■	■		■	■	
	4-110-C5-0	■	■		■	■		■	■		■	■	

Pipettenspitzen Kompatibilität  
AHN Mehrkanal-Pipetten



Informationen entsprechen dem Stand bei Drucklegung. Stand: 29.05.2017



Uthleber Weg 14  
D-99734 Nordhausen  
Germany

Tel.: +49(0)3631/46594-04  
Fax: +49(0)3631/46594-10

E-Mail: [info@ahn-bio.de](mailto:info@ahn-bio.de)  
Web: [www.ahn-bio.de](http://www.ahn-bio.de)