

# **COMPATIBILITY OF**

AHN myTip FT Standard Filter Tips

# AHN myTip FT Standard Filter Tips Compatibility with AHN Single Channel Pipettes

**Filter Tips Compatibility**  
AHN Single Channel Pipettes

Pipettor		AHN																													
Brand		pipet4u <sup>®</sup> pro												pipet4u <sup>®</sup> Revolution																	
Model		variable/single channel						fix/single channel						variable/ electronic single channel																	
Volume Range / Cat. No.		0.1-2.5 µl / 8-100-20-9	0.5-10 µl / 8-101-20-9	2-20 µl / 8-102-20-9	5-50 µl / 8-103-20-9	10-100 µl / 8-104-20-9	20-200 µl / 8-105-20-9	50-200 µl / 8-105-21-9	100-1000 µl / 8-106-20-9	200-1000 µl / 8-106-21-9	500-5000 µl / 8-107-20-9	2000-10000 µl / 8-108-20-9	5 µl / 8-001-20-9	10 µl / 8-002-20-9	20 µl / 8-003-20-9	25 µl / 8-004-20-9	50 µl / 8-007-20-9	100 µl / 8-008-20-9	200 µl / 8-009-20-9	250 µl / 8-010-20-9	500 µl / 8-011-20-9	1000 µl / 8-012-20-9	2000 µl / 8-013-20-9	2500 µl / 8-014-20-9	5000 µl / 8-015-20-9	10000 µl / 8-016-20-9	0.2-10 µl / 8-500-20-0	5-100 µl / 8-501-20-0	10-200 µl / 8-502-20-0	50-1000 µl / 8-503-20-0	100-5000 µl / 8-505-20-0
AHN myTip FT Filter Tips																															
Volume	Cat.No.																														
0.1-10 µl	2-001-70-0 2-001-C5-0	■	■	■																						■					
0.1-10 µl XL	2-061-70-0 2-061-C5-0	■	■	■																						■					
0.5-10/20 µl	2-051-70-0 2-051-C5-0		■	■																											
1-20/30 µl	2-116-70-0 2-116-C5-0				■	■	■	■				■	■	■	■												■	■			
1-50 µl	2-111-70-0 2-111-C5-0				■	■	■	■				■	■	■	■												■	■			
1-100 µl	2-112-70-0 2-112-C5-0				■	■	■	■				■	■	■	■												■	■			
1-200 µl	2-113-70-0 2-113-C5-0					■	■	■								■	■										■	■			
5-300 µl	2-115-70-0 2-115-C5-0					■	■	■								■	■														
100-1000 µl short	2-209-70-2 2-209-C5-2								■	■									■	■	■									■	
100-1000 µl long	2-201-70-0 2-201-C5-0								■	■									■	■	■									■	

Information correct at time of printing. Date: 29 May 17



metal-free

autoclavable  
121° C | 15 min



# AHN myTip FT Standard Filter Tips Compatibility with AHN Multi Channel Pipettes

Pipettor		AHN											
Brand		pipet4u <sup>®</sup> pro				pipet4u <sup>®</sup> Revolution							
Model		variable/multi channel				variable/multi channel							
		8 Channel		12 Channel		8 Channel		12 Channel					
Volume Range / Cat. No.		0.5-10 µl / 8-201-20-9	5-50 µl / 8-202-20-9	3.0-300 µl / 8-203-20-9	0.5-10 µl / 8-401-20-9	5-50 µl / 8-402-20-9	30-300 µl / 8-403-20-9	0.2-10 µl / 8-520-20-0	2.5-50 µl / 8-521-20-0	15-300 µl / 8-522-20-0	0.2-10 µl / 8-530-20-0	2.5-50 µl / 8-531-20-0	15-300 µl / 8-532-20-0
AHN myTip FT Filter Tips													
Volume	Cat.No.												
0.1-10 µl	2-001-70-0	■			■				■			■	
	2-001-C5-0	■			■				■			■	
0.1-10 µl XL	2-061-70-0	■			■				■			■	
	2-061-C5-0	■			■				■			■	
0.5-10/20 µl	2-051-70-0	■			■				■			■	
	2-051-C5-0	■			■				■			■	
1-20/30 µl	2-116-70-0	■	■	■	■	■	■		■	■	■	■	
	2-116-C5-0	■	■	■	■	■	■		■	■	■	■	
1-50 µl	2-111-70-0	■	■	■	■	■	■		■	■	■	■	
	2-111-C5-0	■	■	■	■	■	■		■	■	■	■	
1-100 µl	2-112-70-0	■	■	■	■	■	■		■	■	■	■	
	2-112-C5-0	■	■	■	■	■	■		■	■	■	■	
1-200 µl	2-113-70-0	■	■	■	■	■	■		■	■	■	■	
	2-113-C5-0	■	■	■	■	■	■		■	■	■	■	
5-300 µl	2-115-70-0	■	■	■	■	■	■		■	■	■	■	
	2-115-C5-0	■	■	■	■	■	■		■	■	■	■	

**Filter Tips Compatibility**  
AHN Multi Channel Pipettes





Uthleber Weg 14  
D-99734 Nordhausen  
Germany

Tel.: +49(0)3631/46594-04  
Fax: +49(0)3631/46594-10

E-Mail: [info@ahn-bio.de](mailto:info@ahn-bio.de)  
Web: [www.ahn-bio.de](http://www.ahn-bio.de)